


- Cloning mouse pancreatic α -amylase cDNA into a mammalian expression plasmid

cDNA cloning: Experimental steps

- Isolation of total RNA (Lab 2)
- Reverse transcription (Lab 2)
 - oligo-dT primer anneals to polyA tail of mRNA
 - Reverse transcriptase will synthesize first-strand cDNA
- PCR amplification using first-strand cDNA as a template (Lab 3)
- Ligation of PCR amplicon (Amy2 coding sequence) into expression plasmid / Bacterial Transformation to select recombinant plasmids (Lab 4)


Total RNA Composition



1. RNA isolation
2. Quality control
3. Reverse transcription (RT)

RNAases present on skin etc
degrade RNA!


- Clean benches with 70% ETOH
- Wear gloves
- Use filter tips
- Keep RNA on ice



Total RNA Isolation

- **Trizol extraction** (prepared by us)
 - Guanidinium thiocyanate => Denatures proteins, including RNases, and preserves RNA integrity.
 - Phenol-chloroform extraction => RNA separates into the aqueous phase.
 - [Video](#)
- Purification on **RNeasy spin column** (prepared by you)
- You will isolate RNA from mouse tissue

Trizol Extraction



aqueous phase: RNA

interphase: DNA

organic phase: proteins, lipids



RNA isolation: RNeasy spin column

Use aqueous phase Trizol extract
Selective binding of the RNA to a
silica-based membrane.

Purification of RNAs > 200
nucleotides.

Safety: wear goggles



RNA Purity

NanoDrop Spectrophotometer

Assess purity of total RNA

Compare to control RNA

- Concentration: how much?
- Purity: how clean?
- Integrity: how intact/ how much degraded?

- Nucleic acids have an absorption maximum at at **260 nm**
- The Beer-Lambert law **$A=El$**
- Linear change in **absorbance** with **concentration**
- Spectrophotometric conversion for RNA:
1.0 A₂₆₀ unit = 40 µg/ml (1 cm path)

Chemical Purity of Nucleic Acids is Assessed by Absorbance Ratios


Absorption at 280 nm may be caused by proteins (aromatic amino acids) or phenol

- A260/A280 DNA 1.7 - 1.8
- A260/A280 RNA 1.9 - 2.2


Absorption at 230 nm can be caused by guanidinium thiocyanate, other organic compounds or proteins.

- A260/A230 RNA and DNA > 1.9

- Small samples: 0.5 μl - 2.0 μl
- No need for cuvettes or capillaries.
- With the arm open, a sample is pipetted directly onto the pedestal.
- After the arm is closed, a sample column is formed.
- The pedestal then moves to automatically adjust for an optimal path length (1 mm).
- When the measurement is complete, the surfaces are simply wiped with a lint-free lab wipe before going on to the next sample.



Typical NanoDrop Spectral Profile




← RNA sample

← A260/A230

← A280/A230

← Concentration



RNA integrity

Agarose gel electrophoresis


Assess integrity of ribosomal RNA
bands of purified RNA

Compare to control RNA


- How to run an agarose gel:
 - Wear gloves (GelRed dye)
 - Plug cables correctly (+/-)
 - Turn power ON after loading your sample
 - Watch samples migrate in the correct direction
 - Turn power OFF before taking out the gel



- The gel contains a fluorescent nucleic acid dye (GelRed)
- When illuminated with UV light the nucleic acids are visible
- Save the image with group number



Reverse Transcription (RT) Reaction



- Template: total RNA from **mouse pancreas** (provided by us)
- Oligo-dT primer
- Reverse transcriptase: derived from Moloney Murine Leukemia Virus
- Deoxyribonucleoside triphosphates (dNTPs)

- **Sample without reverse transcriptase (-RT)**

This control assesses the amount of DNA contamination present in an RNA preparation. Important control during PCR amplification (Lab 3)