

CH314 – Structural Analysis Part III : X-ray Tools

Exercise 2

1) Identify the lattice planes



2) Which diffraction patter do you attribute to a) single crystal b) powder c) nanocrystalline powder



3) You perform a powder diffraction experiment on your colloidal gold nanoparticles. From the data, figure out if your particles are spherical or if they have a preferred axis. Consider the Debye-Scherrer formalism.

Journal of Colloid and Interface Science 338 (2009) 105-110



Contents lists available at ScienceDirect

Journal of Colloid and Interface Science





Use of wide-angle X-ray diffraction to measure shape and size of dispersed colloidal particles

S. Junaid S. Qazi a.*, Adrian R. Rennie a, Jeremy K. Cockcroft b, Martin Vickers b

The diffraction experiments were carried out with a Stoe STADI-P X-ray powder diffractometer at University College London with a Ge(111) monochromator before the sample. A Cu Kα₁ X-ray beam of wavelength 0.154056 nm was used. The samples were measured in 0.35 mm diameter glass capillaries at room temperature. Capillary diameters were checked to be equal with a micrometer. The diffraction patterns were collected using a position-sensitive

Lince the particle diameter for various lattice planes is similar the particle appear to be whereal or maybe elliptical.

Table 2 Results from the pseudo-Voigt fit to the Au particles. 2θ is the peak position for the Cu $K\alpha_1$ radiation wavelength (0.154056 nm).

0.114

Instrument

resolution (°)

Example	calculation
V 355	0
0 = 38.3	21

2	. 0 =	= 38,21		31
				42
C=	0 =	19.1050=	12.333 md	

hkl

111

200 44.34 0.822 0.115 220 64.69 0.810 0.123 77.63 0.905 0.135 11 115.34 1.194 0.119

Peak width

'β' (°)

0.580

Bragg angle

2θ (°)

38.21

Make sure to convert to radians!

Corrected width = B = 0.47° = 0.0082 rad

Diameter =
$$\frac{\text{k} \cdot \text{h}}{\text{width} \cdot \cos \theta} = \frac{0.9 \cdot 0.154056 \text{ nm}}{0.0082 \cdot \cos (0.333)} = 17.9 \text{ nm}$$

³ Materials Physics, Uppsala University, Angströmlaboratoriet, Box 530, 75121 Uppsala, Sweden
⁶ Applied Crystallography Group, Department of Chemistry, University College London, 20 Gordon Street, London WC1H 0AJ, UK