


# CH314 – Structural Analysis

## Part III : X-ray Tools

### Exercise 3

1) Lets pick up the problem from last week again. The data below shows diffraction patterns from a solution of colloidal nanoparticles taken with an x-ray tube using the Cu  $1K\alpha$  line with a wavelength of 1.54 Angstrom




**Fig. 1.** X-ray diffraction pattern from  $Ni(OH)_2$  dispersions ( $a = 3.1273(2) \text{ \AA}$ ,  $c = 4.610(3) \text{ \AA}$ ) with diffraction peaks labeled with  $hkl$  indices. The 001 peak is

**Table 1**

Results from the pseudo-Voigt function fit to the  $Ni(OH)_2$  sample.  $2\theta$  is the peak position for the Cu  $K\alpha_1$  radiation wavelength (0.154056 nm).

| $hkl$ | Bragg angle $2\theta$ ( $^\circ$ ) | Peak width ' $\beta$ ' ( $^\circ$ ) | Instrument resolution ( $^\circ$ ) | Corrected $\beta$ ( $^\circ$ ) | Diameter (nm) |
|-------|------------------------------------|-------------------------------------|------------------------------------|--------------------------------|---------------|
| 001   | 19.25                              | 1.140                               | 0.121                              | 1.02                           |               |
| 100   | 33.06                              | 0.197                               | 0.115                              | 0.08                           |               |
| 101   | 38.52                              | 0.593                               | 0.114                              | 0.48                           |               |
| 102   | 52.06                              | 0.876                               | 0.117                              | 0.76                           |               |
| 110   | 59.05                              | 0.248                               | 0.120                              | 0.13                           |               |
| 111   | 62.67                              | 0.550                               | 0.122                              | 0.43                           |               |



**Fig. 2.** X-ray diffraction patterns from dispersed Au particles ( $a = 4.067(3) \text{ \AA}$ ) with diffraction peaks are labeled with  $hkl$  indices.


**Table 2**

Results from the pseudo-Voigt fit to the Au particles.  $2\theta$  is the peak position for the Cu  $K\alpha_1$  radiation wavelength (0.154056 nm).

| $hkl$ | Bragg angle $2\theta$ ( $^\circ$ ) | Peak width ' $\beta$ ' ( $^\circ$ ) | Instrument resolution ( $^\circ$ ) | Corrected $\beta$ ( $^\circ$ ) | Diameter (nm) |
|-------|------------------------------------|-------------------------------------|------------------------------------|--------------------------------|---------------|
| 111   | 38.21                              | 0.580                               | 0.114                              | 0.47                           |               |
| 200   | 44.34                              | 0.822                               | 0.115                              | 0.71                           |               |
| 220   | 64.69                              | 0.810                               | 0.123                              | 0.68                           |               |
| 311   | 77.63                              | 0.905                               | 0.135                              | 0.77                           |               |
| 420   | 115.34                             | 1.194                               | 0.119                              | 1.07                           |               |

- What can you tell about the  $Ni(OH)_2$  and the Au nanoparticles by just looking at the data?
- Calculate the nanoparticle dimensions for the 001 and 100 reflections for the  $Ni(OH)_2$  and 111, 200 reflections for the Au particles.
- Sketch the dimensionality of each of the particles.
- Which lattice do the particles exhibit? Work on the Au nanoparticles first and guess for the  $Ni(OH)_2$  particles.

2) The figure below shows a x-ray diffraction experiment in an Ewald sphere construction. The inverse lattice vector  $\mathbf{a}$  has a length of  $4\text{nm}^{-1}$ .



- Which points will yield a diffraction peak?
- Which Bragg peaks will be observable on the detector?
- What do you need to do to observe the (01) peak?
- What is the photon energy in this experiment? Remember that 1.24 keV corresponds to a wavelength of 1 nm.