



Jigsaw 4B

2D NMR: NOESY and EXSY

1. The 2D NOESY (Nuclear Overhauser Effect Spectroscopy) experiment gives information on through-space interactions, which can be used to identify protons which are close in space. Below, you have the aromatic region of the ^1H 2D NOESY spectrum of 12,14-di t -butylbenzo[*g*]chrysene. Assign the different aromatic ^1H peaks.



2. [Week 3 Slides 71–80] Consider the EXSY spectrum of a symmetrical 2 spin (A and B) system. The intensity of the four peaks are given by the following equations:

$$I_{AA}(\tau_m) = \frac{1}{2}[1 + \exp(-2k\tau_m)]\exp(-\tau_m/T_1)M_{A0}$$

$$I_{BB}(\tau_m) = \frac{1}{2}[1 + \exp(-2k\tau_m)]\exp(-\tau_m/T_1)M_{B0}$$

$$I_{AB}(\tau_m) = \frac{1}{2}[1 - \exp(-2k\tau_m)]\exp(-\tau_m/T_1)M_{B0}$$

$$I_{BA}(\tau_m) = \frac{1}{2}[1 - \exp(-2k\tau_m)]\exp(-\tau_m/T_1)M_{A0}$$

- a. Draw the 2D EXSY spectrum, taking into account the relative intensity of diagonal and cross peaks, at the following values of $k\tau_m$:

i. $k\tau_m = 0$

ii. $k\tau_m \ll 1$

iii. $k\tau_m \approx 2$

iv. $k\tau_m \gg 2$

- b. For each plot in (a), explain what is happening to the peak intensity and why.

- c. 2D exchange spectroscopy is used to determine exchange in which motion regime? Explain why.