

Protein Expression in BL21(DE3) Protocol

1. Resuspend a single colony from the transformation plate in 5 mL of LB media containing 50 $\mu\text{g}/\text{mL}$ ampicillin (AMP) in a Falcon tube.
2. Incubate at 37 °C with shaking at 260 rpm overnight.
3. Inoculate 200 mL of fresh LB media containing 50 $\mu\text{g}/\text{mL}$ AMP with 2 mL of the overnight culture in a 500 mL Erlenmeyer flask.
4. Incubate at 37 °C with shaking at 260 rpm until the OD600 reaches 0.5–0.6 (about 2 hours).
5. Induce protein expression by adding 200 μL of 100 μM Isopropyl β -D-1-thiogalactopyranoside (IPTG) and incubating for an additional 3 hours at 37 °C with shaking.
6. Take 100 μL of the sample for SDS-PAGE.
7. Harvest the cells by centrifugation at 4,000 x g for 10 minutes at 4 °C.
8. For SDS-PAGE:
 - Pellet the 100 μL sample by centrifugation (1–2 min at 13,000 rpm).
 - Resuspend the pellet in 10 μL 2 \times loading buffer + 10 μL MQ water.
 - Incubate at 95°C for 5 min.
 - Briefly spin down and load 8 μL onto the gel.