



The eucentric position

The eucentric position is the horizontal center of the objective lens. The sample must be set to this position. To do this the entire sample holder is raised or lowered.

TEMs are set up so that magnification, camera length, and correct focus are set to this reference position. When the sample height (in the Z-direction) is set at the eucentric position, the sample can be tilted around its axis without the image of the sample moving across the projection screen.

The eucentric height can be set by adjusting the objective lens current to a specific known setting for a specific voltage (e.g. 80kV, 100kV or 200kV), and then raising or lowering the Z-height until the image is in focus. The "wobbler" button can be used as a focus aid where you bring the vibrating split image together into one static image.