

## Subject 6: Electrical and optical characterization


### 1. Objective

During this laboratory, you will learn how to use the four-point probe device for electrical characterization of aluminum thin films and an optical microscope to measure the dimensions of patterns and alignment errors.

### 2. Optical characterization

#### a. Alignment errors

This procedure is used to determine alignment capability using Vernier scales. They are patterns designed to resolve alignment errors for a given process. Figure 1 shows a vernier design. The blue structure has lines that are  $4\text{ }\mu\text{m}$  wide and spaces that are  $3.8\text{ }\mu\text{m}$  wide (resulting in a pitch of  $7.8\text{ }\mu\text{m}$ ). The red structure has both lines and spaces which are  $4\text{ }\mu\text{m}$  wide (the corresponding pitch is  $8\text{ }\mu\text{m}$ ). The difference between the two structures pitch,  $0.2\text{ }\mu\text{m}$ , is the resolution of the vernier. A perfect alignment is obtained when the central lines (long lines) are well aligned. In case a misalignment exists, find where the first red and blue lines are aligned and count its position from the long line. The misalignment is the position number times the offset,  $0.2\text{ }\mu\text{m}$ . A second vernier with a resolution of  $1\text{ }\mu\text{m}$  is added in case where error alignment is very large.



The pattern of the vernier in your mask with features size are as follow:







Figure 1: layout of the Vernier, in red is level one and in blue is the level two. The two levels are shown superimposed.



Example of 0 misalignment.



Example of  $-1\ \mu\text{m}$  misalignment (fifth position times  $0.2\ \mu\text{m}$ ).



Example of  $1\ \mu\text{m}$  misalignment (fifth position times  $0.2\ \mu\text{m}$ )

### b. Laboratory Procedure

Under the optical microscope (Keyence digital microscope), observe your wafer in which you have already aligned the second-level photoresist on the patterned first-level wafer (in the cleanroom session).

- Each student characterizes a vernier located at (x,y) position.
- Inspect the vernier on the wafer and take a picture.
- Record misalignment errors for x and y positions at different locations on the wafer:

Student name				
Vernier coordinates (x,y)	(0,28)	(28,0)	(-28,0)	(0,-28)
$\Delta x$ ( unit )				
$\Delta y$ ( unit )				

- Observe under optical microscope, your wafer with aluminum patterns obtained by the etching process and the one by lift-off process.
- Take pictures of the same pattern with the same magnification.
- Measure their dimensions.
- Inspect and compare the difference between the patterns obtained by the two methods.

### 3. Electrical characterization

#### a. Transmission line method (TLM)

Resistors are tested and characterized with Karl Suss PM6 Probe station and the Keysight 34420A Nano\_Volt/Micro\_Ohm Meter. The devices to be tested are aluminum resistors obtained by **chemical etching** and **lift-off** methods. The layout of the resistors is shown in figure 2.





Figure 2. layout of resistor patterns

These structures will be used to plot the curve  $R$  versus  $L$ , where  $R$  is the electrical resistance and  $L$  is the length of the resistor segment.

#### b. Laboratory Procedure

- Place the sample on the wafer stage




- Place the probe tips on the contact areas as follows:
  - Two probe tips for current injection, these probes are labeled  $I_1$  and  $I_2$ . Under the optical microscope, use the micropositioning stage to place each tip on the corresponding pad on your sample with the same label of the probe tip.
  - Two probe tips for voltage measurement with labels  $V_1$  and  $V_2$ . Under the optical microscope, place each tip on the corresponding pad on your sample.
  - Select the 4W mode on the Keysight 34420A Nano\_Volt/Micro\_Ohm Meter, the resistance value is directly displayed.
  - complete the following table:

contact	V1-V2	V1-V3	V2-V3	V2-V4	V3-V4	V4-V5	W (unit)	Student name
L(μm)	400	750	300	550	200	100		
R(unit)								
R(unit)								
R(unit)								
R(unit)								

- Plot R versus L and calculate the slope of the curve.

c. *Sheet resistance measurement*

The following Structures will be used to measure the sheet resistance.



d. *Laboratory Procedure*

- Place the probe tips on the contact areas as follows:
  - for current injection: pads 1 and 2
  - for voltage measurement: pads 3 and 4.
  - The resistance value displayed is used to calculate the square resistance by the following formula:

$$R_s = \frac{\pi}{\ln 2} \times \frac{V_{34}}{I_{12}}$$

- Complete the following table :

structure	A	B	C	D	E	F	G	H
Rs (unit)								


- Compare the value of Rs for the different structures.
- From Rs and the slope of the curve R versus L, calculate the width of the measured structure.
- Compare this value to the one obtained by optical measurement.
- Compare the measured value of W to the one registered on the photomask.

- Knowing the thickness of the aluminum film (measure done in the cleanroom session), determine the resistivity of the film?

EPFL

## Résistance électrique

### Piste de section rectangulaire



$$R = \rho \frac{L}{W \cdot T}$$

$$L = W \rightarrow R = \frac{\rho}{T} = R_{\square}$$

$$R = R_{\square} \frac{L}{W}$$


L (longueur) et W (largeur): caractérisent la géométrie du masque

$R_{\square}$  Résistance par carrée : caractérise la technologie

$\rho$  Résistivité du film


T Epaisseur du film


### Barre de section rectangulaire



$$R = R_{\square} \frac{L}{W}$$




$$pente = \frac{R_{\square}}{W}$$





$$R_{\blacksquare} = R_{sh} = \frac{\pi}{\ln(2)} \frac{V_{34}}{I_{12}}$$

$$V_{45} = \frac{R_{sh} L I_{26}}{W}$$



$$R_{\blacksquare} = R_{sh} = \frac{\pi}{\ln(2)} \frac{V_{34}}{I_{12}}$$

$$V_{45} = \frac{R_{sh} L I_{26}}{W}$$

