




Towards large (bio)molecules




- Low concentrations
- Many vibrations

Photo fragmentation ion spectroscopy



Background free technique




UV photodissociation spectroscopy



$T = 6 \text{ K}$




Gas-phase but solution-like






- UVPD produces charged fragments

UV photodissociation spectroscopy of cold ions




$T = 6 \text{ K}$





- Cooling enables vibrational resolution
- **UVPD** produces charged fragments



IR-UV photodissociation spectroscopy of cold ions




$T = 6 \text{ K}$





- Cooling enables vibrational resolution
- UVPD produces charged fragments





***IR-UV** photodissociation spectroscopy of cold ions*



T = 6 K




Conformer selective IR spectroscopy





Depletion

Three distinct conformers



Gramicidin-S: validating theory




Natural antibiotic:



10 amino acids, 176 atoms.

Experiment and accurate theory are limited by size of molecules

Conformer selective UV spectroscopy




**Integral
fragmentation is
conserved**



Absolute *IR* absorption cross-section

Because $h\nu_{\text{UV}} \gg D_0 \gg h\nu_{\text{IR}}$: $\beta \simeq \text{const.}$



Fix *UV*, Scan *IR*

UV absorption:

$$I_A(\nu) = n_0 \sigma_A^0 \Phi^{UV} + n_1 \sigma_A^1 \Phi^{UV}$$

Population:

$$n_0(\nu) = n_p \exp(-\sigma^{IR}(\nu) \cdot \Phi^{IR})$$

$$n_1(\nu) = n_p [1 - \exp(-\sigma^{IR}(\nu) \cdot \Phi^{IR})]$$

$$n_0 + n_1 = n_p$$

UV absorption probabilities:

$$\Phi^{IR} = 0 \Rightarrow I_A^0 = n_p \sigma_A^0 \Phi^{UV}$$

$$\Phi^{IR} \sigma^{IR} \rightarrow \infty \Rightarrow I_A^\infty = n_p \sigma_A^1 \Phi^{UV}$$

$$\sigma_A^{IR}(\nu) = \ln\left(\frac{I_A^0 - I_A^\infty}{I_A(\nu) - I_A^\infty}\right) \frac{1}{\Phi^{IR}}$$

Analytical application of photofragmentation ion spectroscopy

For what?

Distinguish isomeric biomolecules

Why?

Biological properties depend on isomers; many biomolecules are highly abundant in isomers.

How?

By comparison of their spectra and/or combined spectroscopic/mass-spectrometric fingerprints: library-based approach.

What should be done?


- Measure and store the spectra of all known isomers (library of standards)
- Measure under identical conditions a spectrum of an unknown isomer or a mixture of them
- Make linear numerical fit of the analytical spectrum by the sum of the library standards.

$$A(\lambda) \simeq \alpha_1 \cdot S_1(\lambda) + \alpha_2 \cdot S_2(\lambda) + \dots + \alpha_n \cdot S_n(\lambda);$$

$$\alpha_1 + \alpha_2 + \alpha_3 + \dots + \alpha_n = 1.$$

- The coefficients of the best fit will give the relative concentrations of the isomers.


Identification of solvent-tagged biomolecules by **IR** spectroscopy





Produce ions retaining a few solvent molecules (H-bonds);

Measure IR photodissociation spectra by detecting loss of solvent molecules

Apply library-based approach





UV photodissociation spectroscopy of cold ions





Boyarkin et al., JACS (2006)

- Cooling enables vibrational resolution
- **UVPD** produces charged fragments
- Orbitrap **MS** detects all fragments at once





Colors for cold masses: 2D UV-MS fingerprinting

$T = 6\text{ K}$



2D UV-MS fingerprint


- Extremely sensitive to structural details
- Reflects fundamental property ($\Delta E, m$) =>
- Highly reproducible, low sensitive to exp conditions



2D UV-*MS* Library-based quantification of isomers

Unknown
Solution mixture

Library of standards (4 isomers):






α_i – relative concentrations of mixed in solution molecules;

Identification of components in a mixture of isobaric peptides


Isomeric phosphopeptides:

TSAAATSpY

TSAAApTSY



UV photofragment spectroscopy: enkephalin



YAGFL


([Ala²]-Leu-enkephalin

Low potent pain relief drug


[D-Ala², D-Leu⁵]-enkephalin
(DADLE)

High potent pain relief drug

The ions are cooled
to T= 10 K



UV photofragment spectroscopy: enkephalin



YAGFL


([Ala²]-Leu-enkephalin

Low potent pain relief drug

[D-Ala², D-Leu⁵]-enkephalin
(DADLE)


High potent pain relief drug

The ions are cooled
to T= 10 K




Quantification of isomeric ephedrines

- Central nervous system stimulant (doping) prohibited in sports by WADA:






Ephedrine:

$$11 \pm 0.5 \text{ } \mu\text{g/mL}$$


Pseudoephedrine:

$$170 \pm 0.5 \text{ } \mu\text{g/mL}$$


- Analog of substituted β -phenylethylamines (forbidden psychoactive drugs)

R=CH₃/OH/H

Isomeric glycans


Monosaccharide pyranose :

$2^5 = 32$ isomers



Tetra-saccharide:

35 560 isomers....




Types of isomers:

- 1) α/β – anomers
- 2) epimers (Gal, Glc)
- 3) connectivity isomers
- 4) configuration
- 5) structural
- 6) D/L enantiomers




Breast milk:
>200 glycans

Nestlé formula:
< 10 glycans



D-glucose




Galactose




Both are sweet, but only **D**
(natural) gives energy;
L cannot be metabolized

Do NOT absorb in UV !!!

UV non-absorbing biomolecules: isomeric monosaccharaide



Strong non-covalent interactions in glycan...Tyr complexes



GalNac

GlcNac

Strong interactions


✓ $\text{HO}_{\text{Tyr}} \cdots \text{HO}_{\text{Carb}}$ –
Lowers energy of the
oxygen n -orbital:
blue shift

✓ $\text{OH}_{\text{Tyr}} \cdots \text{OH}_{\text{Carb}}$
Inductive polarization
of the π -system:
red shift


Weaker interactions...

Libraries of isomeric mono- and disaccharides

TrmH⁺--R



TyrH⁺--R




UV photofragment spectroscopy for identification of oligosaccharides

<i>Real.</i>	<i>Calc.</i>	<i>Real.</i>	<i>Calc.</i>
<i>Hepta-II</i>		<i>Hepta-I</i>	
66,7%	71%	33,3%	29%


2:1 and fit

Heptasaccharides
(174 atoms)



Identification of what is not yet available

Enantiomers (cannot be distinguish by absorption)



Predicting 2D UV-MS identity of a **nonexistent** enantiomer

Laser Chemistry

□

- Selective heating of reactants only;
- Selective heating of one sort of molecules only
(chemical species, isomers, isotopes);
- Selective excitation of vibrations only;
- Collisional free dissociation;
- Selective excitation of a specific chemical bond

Isomerisation

$$K = \exp(-\Delta G / RT)$$



1.2-dichloroethene

cis-

trans-

hexafluorocyclobutene

hexafluoro-1,3-butadiene



Isomerisation

7-dehydrocholesterol

7-dehydrocholesterol



248 nm

pre-vitamin D3




248 nm

90% yield


vitamin D3

2% yield


Side reactions are suppressed
because of internal heating only.



(I)



(II)




(III)

Examples of laser induced reactions



- Polymerization through a pulsed creation of radicals:

Chain length is controlled by change of dark time;



- Laser pyrolysis of SiH_4 for Si deposition:

Deposition temperature is controlled independently of dissociation;



No surface heating




Those are collisional single photon processes

Infrared Multiple Photon Excitation and Dissociation


Anharmonicity compensation by:

- **Intermode anharmonic couplings; (IVR)**
- Degeneracy lifting;
- Two-three photon transitions;
- Power broadening;
- Increase of rotational energy (P $\rightarrow$ Q $\rightarrow$ R route);

Vibrational anharmonicity




Infrared Multiple Photon Excitation and Dissociation





Regardless of which bond is excited,
the weakest bond always breaks first!

IR laser chemistry



T = 425 K,
t = 20 h.



Room t° ,
t = 5 min.




50%

Molecular Laser Isotope Separation

Isotopes in nature

Oxygen



O-13	0.0089 s
O-14	70.599 s
O-15	122.24 s
O-16 stable, 99.762%	
O-17 stable, 0.038%	
O-18 stable, 0.200%	
O-19	26.91 s
O-20	13.57 s
O-21	3.4 s

^{13}C natural abundance:

on Earth: ^{13}C - 1.07%, ^{12}C - 98.9% (on some stars ^{13}C - 22%).


Nonlaser isotope separation techniques:

GAS DIFFUSION

$$\bar{v} = \sqrt{8RT/\pi M}$$

$$K(^{235}\text{U}) = 4 \cdot 10^{-3}$$


- Requires many cycles,
- very high energy consumption.




GAS CENTRIFUGE

$$n(r) = n(r_0) \cdot \exp(-m(r_0 - r)^2 \omega^2 / 2kT)$$


Lighter gas has higher concentration closer to the axis



Cryogenic distillation



ISOTECH PLANT



Principles of laser isotope separation methods

- Selective electronic/vibrational excitation of targeted isotopic species,
(Based on isotopic shift)
- Conversion of excited species to chemically/physically new species without loss of selectivity,
(Ions, new sort of molecules, clusters, etc)
- Physical separation of the new enriched species from parent ones
(temperature distillation, electromagnetic separation of ions and neutrals, etc).

Methods

1. Atomic LIS
2. Molecular



Photoionization approach to LIS

Basis: shift of electronic energy levels due to different number of neutrons




Different frequencies of atomic spectral transitions


*Doppler free U atom
absorption spectrum.*



Photoionization approach (AVLIS)







Chain of dye lasers at AVLIS pilot plant




AVLIS facilities at LLL




Molecular approach for ^{13}C isotope separation (MLIS)




10-15 kg/year, 30-50% of ^{13}C .




Overtone preexcitation-IRMPD approach



CF_3H^*




3_1 band




3_1 CH stretch-bend vibration in CF_3H_{43}

Overtone preexcitation-IRMPD approach





Collisional effect in OP-IRMPD of CF₃H at $\nu_{\text{CH}}=3$:

pressure dependence



Up to 16-fold increase
of selectivity by collisions



Traditional vs Laser methods

<i>Traditional technologies</i>	<i>E, eV/atom</i>	<i>Laser approaches</i>	<i>E, eV/atom</i>
Gas diffusion, ^{235}U	$3 \cdot 10^6$	Photoionization of atoms, ^{235}U (AVLIS)	10^3
Gas centrifuge, ^{235}U	$4.5 \cdot 10^5$	IR MPD of CF_2HCl , carbon-13 to 95%	10^3
Distillation, ^{10}B	10^3	Overtone preexcitation - IRMPD of CF_3H (developed in EPFL) ^{13}C to 98%	10^2