## **The Procedure**

Patient (secured in a rescue bag) and doctor are lifted, patient with feet toward the front of the helicopter.

The patient and doctor are winched to the maximum height.

The winch is retract to its flying position


The winch operator extends the AeroGlide

While loosening the winch, the winch operator assist the doctor in rotating the patient 90°, his feet enter the helicopter first.

A safety line is attached to the patient and the doctor. The winch is unhooked and raised again to its maximum height.

Using the platform's rollers, the patient is slide onto the stretcher.

Once the patient is safely on the stretcher, the platform is retracted back to its stowed position.



Integrated Aerial Patient Transfer System for Helicopter Rescue Operations





Raphael JAVARY Corentin BARUT
Romain LATTION Rémy GRANGE


## **AeroGlide**


Lightweight retractable rescue platform

The AeroGlide offers a simple yet highly effective solution to the complex challenge of mid-air patient transfers. Its streamlined design ensures smooth, safe transfers, with the platform extending enough to position the winch bag's center of gravity securely. The rollers, featuring anti-backward locking, allow controlled sliding of the patient. Easy to mount and dismount, the AeroGlide requires minimal maintenance. Lightweight, retractable, and cost-effective, it combines safety and reliability for rescue operations.



## **Mechanism**





**AeroGlide in Flight Position** 

**AeroGlide in Operational Position** 

Minimum installation and modification

Light weight, 18.5kg

Robust element
Wistand 200kg (securtiy factor 3)

Space-efficient Ground-level air 0.196m<sup>2</sup>

**Easy to manouver** 

**Easy to install** 

The simple, reliable, and cost-effective solution that ensures safe and efficient mid-air patient transfers, designed to perform when every second counts.



# Robotic Integration for Safe Extraction

ENHANCE SAFETY,
EFFICIENCY, AND
COMFORT WITH RISE.
DESIGNED SPECIFICALLY
FOR THE REGA H-145 D3



## **GROUP 2**

Max A. Aebi Raphaelle Paracuellos Jan Röthlisberger Thomas Roulet

Assist:
Konstantin Jacoby
Professor:
Jürg Alexander Schiffmann

## SEAMLESS PATIENT TRANSFERS IN FLIGHT

- Eliminates the need of landing for patient rescue
- Streamlines operations, decreasing rescue times.
- Minimal crew retraining required by using existing tools.
- Built to endure tough environments
- Ensures smooth transitions, reducing patient discomfort



## ENGINEERED FOR EXCELLENCE


- Secure doctor to hoist inside the helicopter and engage robotic arm drive to move doctor outside of the helicopter
- Descend and secure patient
- Lift patient
- Engage robotic arm drive to move patient back inside
- Secure patient inside the helicopter






### **Technical data**

- Lift up to 270kg
- Available for H145
- Transfer time of 20s
- Designed to be compatible with all currently used systems
- Patient swinging to angles of more than 45° can be handled by system





## MISSION SUCCESS, EVERY TIME

Get in touch today to learn more about how we can improve your patient transport operations. Project Manager Max A. Aebi is available to answer your questions under max.aebi@epfl.ch photo credit: <a href="mailto:www.rega.ch">www.rega.ch</a>

Scan this QR for animations of the transfer procedure



#### Insertion in the H150



The stretcher (fixed on the roller in (2)) then rotates and slides inside of the helicopter on the rail, with the copilot door open.


#### Where to fix it?




Using existing fixation points: two on the central line of the ceiling, one next to the side door on the pilot side and another one next to the copilot's side door.

## **Security**

 Safe transfer: a blocking insures the fixation of the rail with the winch attachment



- Doctor and copilot safety: They do not get exposed outside of the helicopter, they handle the system from the inside.
- Minimal training required.
- No power consumption so no effect on the fuel of the helicopter or compromise of other systems.





## Our concept: Sliding rail system

A sliding rail comes out of the helicopter and blocks itself to a winch attachment. The stretcher, fixed on a roller in the winch attachement then slides on the rail and enters the helicopter. The copilot's door opens to facilitate stretcher rotation for seamless entry. Once inside, the stretcher is detached, the rail comes back inside, and the helicopter resumes flight. This system ensures a fast, safe, and efficient mid-air patient transfer.

#### Our system is...



fast: 4 minutes max transfer



CG approved



0 W



170 Kg




3800 CHF

#### **Minimum modifications**

- The system is fixed on built-in fixation points of the helicopter requiring no need for additional attachments.
- It takes no ground space as it's attached to the cab ceiling.




#### How does it works?



The rail (1) is fixed on a sliding mechanism (3), it slides outside and connects to a winch attachment (2) which has the stretcher attached inside.

#### The mechanism



The winch attachment (2); box fixed on the hoist cable, contains a rolling mechanism that allows the stretcher to slide along the track. The rail fixes itself to this box.

## **Our Strenghts**

SkySafe™ Medical Arm is cleverly designed to support all rescue scenarios: reducing response time, minimizing risks, and improving access to difficult locations. Its efficient and practical approach enhances safety and reliability in critical patient transfer operations.



#### **SAFETY**

Ensuring the patient's safety during the transfer is our main concern. Our design allows to eliminate the risk of injuries for both patient and staff.



#### **EASE OF USE**

The system is designed to ensure ease of use and intuitive movements throughout the steps for easy transfers.



#### STRUCTURAL MODIFICATIONS

The implementation of our system requires very little structural modification and benefits from using the existing winch, reducing the time required to upgrade helicopters.

# Detailed Characteristics



11 kg



#### **Duration**

6 min



**Max Load** 

150 kg



0.975 m

# SKYSAFE<sup>TM</sup> MEDICAL ARM

GROUP 4
Simon Barbe
Simon de Bazelaire
Manon Oulevey
Audrey Ventura







Rega plays a vital role in delivering emergency medical aid across Switzerland's challenging terrains, relying on the Airbus H145 helicopter. Current rescue procedures, which involve landing at intermediate sites, are time-consuming and complex. The SSMA proposes an Integrated Aerial Patient Transfer System to enable secure, uninterrupted patient transfers directly into the helicopter, saving time, reducing risks, and improving access efficiently.

## Our system

SkySafe™ Medical Arm is a 3degree of freedom mechanism, consisting of a horizontal articulated arm. Reaching the patient outside, the arm allows to transfer the patient safely inside, directly on the medical stretcher. The system is integrated inside helicopter, the above the stretcher.




## **Detailed steps**

#### 1 - Bring up and attach to arm




Retract fully the hoist cable and take the arm's ropes from the rescue aid and attach them to the arm.


#### 2 - Transfer the patient inside



Rotating the arm, transfer the patient carefully inside. Move it until the final position above the stretcher is reached.

#### 3 & 4 - Bring down & secure the patient






Using Petzl descender, cautiously bring down the patient on the stretcher. Then secure the patient on the stretcher with the belt system. Store the arm next to the window.

### **Performance**

## **Usability**

The system has been designed to be easily used by the paramedic while airborne. It is compact and compatible with other onboard systems (external hoist, rescue bag, electrical power, ...)





## Size - Weight

The arm fits inside the cabin while allowing the paramedic to move. It weighs less than 35 kg and can lift a 120 kg patient.

## **Safety**

The system ensures the safety of the patient and the emergency flight physician while transferring them into the helicopter. It allows for the sequence of previous operations in case of any undesired problems.





## **Speed**

This system is quick and does not require many manipulations. The streamlined process makes for a swift and efficient rescue.



**Group 5** 



# Integrated Aerial Patient Transfer System for Helicopter Rescue Operations





## SAFETY

Ensure the well-being of all parties, at all times.



Ensure the possibility to deploy any existing rescue mechanisms.

#### ) EASE OF USE

Streamline the procedure for efficient patient rescue.




## **(**


#### - COMPATIBILITY

Maintain functionality across all onboard systems and equipment.






The curved trajectory of the rail, drawn in black, makes the patient move towards the back of the helicopter, over the medical table. The trajectory is shown here in red.



## **Our System**



#### 1. The Table

The rail on which the table slides is extended. It has two secure positions, front or back.


#### 2. The Arm

The arm moves along the rail on its wheels. It features a small electric hoist

#### 3. The Rail

The rail guides the arm's wheels along a curved trajectory.

- The table is slid into its front position
- The patient is brought up to the cabin's level by the outside hoist
- The winch folds towards the cabin
- The arm slides outside the cabin
- The patient is transferred to the arm's carabiner
- The arm slides back inside, bringing the patient over the table
- The patient is lowered onto the table using the inside hoist
- The patient is secured, and the table slides back for transport







# THE JET ANGEL

AERIAL PATIENT TRANSFER
SYSTEM FOR HELICOPTER
RESCUE OPERATIONS




#### Our team:















#### Total cost:

45'420 CHF

(per system considering the production of 20 systems)





## **₹**5

#### Rescue without landing

Eliminates the need to land to transfer patient inside



#### **Ease of operation**

Intuitive process that requires minimal effort from the staff



#### **Easy installation**

This product requires minimal modifications to the helicopter


## Why it works

#### Component specification

- The **rope** is 5m long and can carry 400 kg.
- The carabiner holds up to 2549 kg.
- The **motor** has a nominal torque of 86.55 Nm and a nominal speed of 0.2 m/s.
- The winch is fixed to two fixation points inside the helicopter that can hold up 200 kg each. (This information was found in the flight manual.)

#### Quality assurance

- We made many calculations:
  - Required length of rope
  - Drum dimensions based on rope
  - Center of mass change of the helicopter
- We then added a security factor to those numbers
- We also **compared** the bought parts, like the motor or the bearings, to find the one with the best capabilities and adapted to our needs.
- We did some Finite Element Methods (FEM) simulations on the drum and fixations to make sure it can hold the expected forces.



## More infos

#### **Advantages**

- User friendly
- Minimal number of components
- Many parts are available for purchase
- Lightweight
- Installation requires no modification of the cabin

#### **Potential optimization points**

- Improving ergonomics
- Refining fixation geometry
- Reducing overall weight
- Lowering the cost of parts
- Optimizing size

—> With a bit more time these optimization would have been implemented.

#### Global data

Weight: 15kgVolume: 0.01m^3

Dimensions: 532 x 119 x 196 mmPrice (estimated): ~10'000CHF

• Time to bring patient inside : 40s to 1min

#### Thank you for your attention!

Do not hesitate to read our report if you you want to find out about our calculations and simulations or to see the assembly steps and more!








# Integrated Aerial Patient Transfer System for Helicopter Rescue Operations

ME-403 Group 7

Tania Rosset Bron Ward Alexis Ruprecht Joshua Levy-Baron

supervised by Anupam Jena




## **Our Solution**

<u>What?</u> To transfer the patient inside the helicopter we want to add a **secondary** smaller winch.

Where ? It is inside the helicopter at the opposite of the door.


<u>How it works?</u> A DC motor connected to a reducer makes a **drum rotate**. A rope is rolled up around this drum and it can be unrolled as wiched.

<u>How to use?</u> Controlled by a **remote controller**. It carries the patient inside with the help of the outside winch.







The patient is safely inside the stretcher (red triangle)




How to use it






Step 3: Someone attaches the carabiner of the inside winch to the patient

<u>Step 2:</u> The doctor attaches one of his two carabiner to the inside of the helicopter and **walks in** 




Step 4: The doctor and staff pushes and rotate the patient inside. He is lifted by both winches.



<u>Step 5:</u> Once the patient is on the stretcher, someone can **detach** the patient from the **outside winch**...
<u>Step 6:</u> ... and then from the **inside winch**.



## Integrated Aerial Patient Transfer System for Helicopter Rescue Operations



#### Overview

To counter REGA's current issues, our team proposes a solution to reduce time loss during aerial rescue missions.

Our design will enable REGA to transfer patients from the ground to the helicopter without landing.

## Secondary winch system

We have installed a secondary winch on the opposite side of the helicopter. This system allows the patient to be pulled from outside to inside the helicopter through the side door. Once the patient is secured, the attachment point is fixed to the ceiling to avoid any difficulties with crew movement in the helicopter.





#### Remote control

To ensure an easy guidance for the attachment point, we propose a synchronous motion control between the two winches. Hand and pedal remotes allow an easy displacement in the x and y axis (winches - stretcher plane) as well as independent control on each winch.

## Railing system

Our system comprises a rail fixed to the floor of the helicopter and sliding blocks mounted on the stretcher,

guaranteeing perfectly guided and safe movement for the patient. Thanks to their flared shape, the rails allow rapid attachment of the stretcher, which is held in the grooves by a spring mechanism. Integrated notches lock the stretcher in the transport position and allow it to be easily unlocked by pressing.



## Patient safety

Patient safety is our priority. The secondary winch and rail system ensure complete guidance of the patient during transfer to the helicopter. The self-locking sliders provide a hands-free experience for the helicopter crew, allowing for quicker care of the patient.

## **Implementation**

The design that we propose uses the support pins already present on the helicopter for the secondary winch holding structure, and keeps the overall helicopter structure intact. Minor modifications are needed on the ceiling, to provide an opening for the secondary cable, the attachment point holder and the redirection pulley.

## **ABOUT US**

Imagine a life-saving mission where **every second counts**, but traditional systems make loading patients mid-flight nearly impossible.

Our company changes this with an innovative solution that safely integrates stretchers into helicopters while airborne.

This breakthrough simplifies rescue operations, ensuring **faster & safer** patient care in the moments that matter most.





## SKYRGENT


"Where security meets safety"


## **OUR OFFER**

### Speed and Efficiency in Emergency Situations

With an ergonomic and intuitive design, patient loading is performed in record time, crucial for rescue operations where every second matters.







#### Optimal Safety for Patients and Crew

Secure locking mechanisms provide quick and reliable fixation

minimizing the risk of movement or falls during critical maneuvers.

#### **Robust and Durable Design**

Built with high-quality, lightweight materials, the system combines durability with efficiency. Its low-maintenance construction guarantees long-term usage while reducing operational downtime.

## **PRODUCT**

### **Specifications:**

**Unitary Price** 

10'000 CHF

Expected lifespan

10 years

Rescue time

20 min

Load capacity

170 kg

Device Total Weight

47kg

#### **Advantages:**

Ti

Time saving



**Money saving** 



Allows longer flights



Benefit after 2 years



## Integrated Aerial Patient Transfer System



IAPTS: a removable device designed for Airbus H145 helicopters. Enables the rescue of individuals via helicopter and their transfer into the cabin during hover flight. Thus eliminating the need for an intermediate landing, saving both time and fuel.

BOOST EFFICIENCY, SAVE ENERGY, MISSION'S READY!

#### **DEVICE WEIGHT:**

24.5 kg

#### **MAXIMUM LOAD:**


270 kg

#### LOW COMPLEXITY:

Low part count, compatible with REGA's winch system

#### **REMOVABLE:**

No major modifications on the helicopter for installation



## Pulley, Hinged opening

Easy opening,

Pull the cable inside

# Voltage: 12 - 72 VDC Speed:

250 rpm - 1250 mm/min

#### **WAYS OF IMPROVEMENT:**

Folding mechanism → save space in cabin
Optimisation of part shape → Weight reduction
Modular components → Adaptability to other model