

Prof. Hatice ALTUG

BioNanoPhotonic Systems Laboratory


École Polytechnique Fédérale de Lausanne, Switzerland



What is a biosensor?




Biosensor is an analytical device that can detect and analyze biochemical analytes



Detected Analytes:

- Small molecules (i.e. drugs)
- Biomolecules (i.e. enzyme, hormones)
- Proteins, antibodies
- Nucleic acids (DNA, RNA, miRNA)
- Extracellular vesicles (exosomes)
- Pathogens (virus, bacteria)
- Cells ...




Test/Device Response:

- Yes/No indicator
- Quantitative detection
- Analyze interactions & kinetics
- Extract chemical/conformational information ...

Biosensors have a wide range of applications



Life Sciences



Pharma




Diagnostics



https://upload.wikimedia.org/wikipedia/commons/6/60/IL1b_Crystal_Structure.png

https://www.contractpharma.com/issues/2017-03-01/view_features/modeling-simulation-for-drug-development-formulation/

http://www.pmlive.com/blogs/digital_intelligence/archive/2016/june/roche_launches_bluetooth_coagulation_tes



Detected Analytes:

- Small molecules (i.e. drugs)
- Biomolecules (i.e. proteins, antibodies)
- Nucleic acids (DNA, RNA, miRNA)
- Extracellular vesicles (exosomes)
- Pathogens (virus, bacteria)
- Cells

Enviromental



Food Safety



Security



https://www.gettyimages.ch/detail/foto/infrared-television-display-lizenzfreies-bild/1150080006?adppopup=true https://www.prescouter.com/2018/01/non-invasive-biosensing-companies-2018/

https://primefeed.in/news/4470654/global-environmental-monitoring-sensors-market-2020-competitive-analysis-advantech-crouzet-analog-devices-inc-alps/

BIOSENSORS in DIAGNOSTICS

- Diagnostics is the FIRST step in medical care
- Follow-up treatments are based on diagnosis

PCR Technology



Advantageous:

- Provides accurate results
- High-throughput

COVID-19 EXPERIENCE!

Limitations:

- Time consuming
- Centralized labs
- Expensive instrumentation
- Trained personnel
- Not available to everyone *resource-poor settings*

Rapid Tests (Lateral Flow Assays)



Advantageous:


- Point-of-care technology
- Low-cost
- Used by anyone
- Provides rapid results

Limitations:

- Not accurate/sensitive
- Not quantitative
- Limited multiplexing capability


Diagnostics - Tomorrow








Global health care
Point-of-care diagnostics
Personalized medicine




Diagnostics - Tomorrow





- ☐ Sensitive
- ☐ Throughput
- Multiplexed
- ☐ Real-time
- ☐ Low-cost
- Rapid
- ☐ Integrated
- Portable
- ☐ Intelligent
- ☐ Autonomous


Global health care
Point-of-care diagnostics
Personalized medicine



Next-generation biosensors will be everywhere!



- ☐ Sensors that can measure everything and everywhere will serve at the nodes of the "Internet-of-body" & "internet-of-things"
- ... connected devices enable smarter systems



MAIN BIOSENSOR PARAMETERS




- ☐ Response:
- ...Yes/No or Quantitative




- ☐ Selectivity
- Affinity based

(surface functionalization)


☐ Reliability (Accuracy)



☐ Limit of Detection (LoD) & Dynamic Range fg/mL, pg/mL, ng/mL, mg/mL..



- Concentration
- ☐ Multiplexing Microarray




☐ Throughput

How many samples per test?

☐ Cost
Reader, cartridge

☐ Time Resolution & Response Time
End point, Real-time (Kinetics) or Continuous





☐ Sample Volume

... pL, nL, uL ...


- ☐ Automation & Sample Handling Microfluidics ..
- ☐ Single-Use (Disposable) or Multiple Use Cartridge? Fixed Reader?
- ☐ Footprint & Integration

Bench-top, Portable, Wearable..


Optical Transduction Mechanisms



Biosensor is an analytical device that can detect and analyze biochemical substances








IR absorption


UV-Vis absorption

scattering



Raman scattering

fluorescence




Bioimaging & Sensing

Refractometric detection

Hyperspectral imaging

WHY NANO-PHOTONICS ?





Photonics provides:


- Immunity to electromagnetic interferences
- Ultra-wide bandwidth ...

But, micron long optical wavelengths introduces size mismatch between light & nanometric samples:

Lower performance


WHY NANO-PHOTONICS ?






Optical nano-resonators


- Focus light into nanoscale volumes
- Enable strong light-matter interaction
- Ideal to achieve high sensitivity, small footprint, multiplexed integration, ...





Localization of E&M field in nanophotonic structures






Size comparison between common biological analytes and the decay length of the evanescent field in various nanophotonic structures

Integration Schemes




Vertical Integration



Planar Integration


Commonly Used Label-Free Detection Principles





Resonance engineering

- High field enhancement
- Spatial overlap between analyte & enhanced near-fields
- Q factor
- Resonance wavelength




Nanophotonic Structures




Control optical response and functionality of nano-resonators with:

- Material choice
 - Metals, 2D materials, dielectrics, ...
- Shape, size & arrangement of nanostructures


Metallic Plasmonics




PNAS 2009



Small 2018





Nature Materials 2012



Nature Comm. 2018


Dielectric Metasurfaces






Science 2018, Science Advance 2019 Nature Photonics 2019 Nature Nanotechnology 2019


2D Quantum Materials





Tunable Mid-IR Sensors Science 2015 Light Science & App 2017

Hybrid System



Combine different materials & systems

Advanced Functional Materials 2020 Nano Letters 2020


BIOS: Research Directions



Point-of-Care Diagnostics



- Low-cost, portable, rapid & multiplexed microarray technologies.
- New device schemes including lens-free microscopy, hyperspectral imaging ...




PNAS 2010 ...

Light Science and Technology 2014 & 2018 ACS Photonics 2015, patent granted ACS Nano 2018 Nature Photonics 2019 Small 2020 Nature Communications 2021 ...

Surface Enhanced Mid-IR Spectroscopy



- Chemical-specific detection
- Conformational sensitive



PNAS 2009 ...

Nature Communications 2013, patent granted Science 2015 & 2018 Nano Letters 2016, 2018 Light Science & Applications 2016 & 2017 Nature Communications 2018 ACS Sensor 2018 ACS Photonics 2018 Advanced Materials 2021 & 2022

Science Advances 2019 & 2023 ...

Live Cell & Single-Cell Technologies



- · Cell screening & profiling, e.g. for immunoengineering applications
- Real-time & label-free secretion analysis from single cells
- Cell-cell interactions & incorporation of biomimetic surfaces











Lab-on-a-Chip 2017 ACS Sensors 2018


Small 2018


Science Advances 2021


Biosensors and Bioelectronics 2022 Nature Biomedical Engineering 2023 ...


11

Infrared spectroscopy provides chemical specific identification










- Chemical specific
 - Label-free





Mid-infrared sensing for biology



non invasive access to molecular vibrations

- signal intrinsic to molecular structure
- Resolves secondary structure from amide I band shape






Challenges of Infrared Spectroscopy



low sensitivity


Beer's law
$$T = \frac{I}{I_0} \sim \exp(-\alpha d)$$

- Prohibits measurements from monolayers with nanometer thickness
- Requires extremely high concentrations



bulky & expensive instrumentation


State-of-the-art is the table-top FTIR



water absorption

Water bands **overlap** with the **vibrational signatures of biomolecules**


 Requires sample drying or D₂0 which are not compatible with living systems (i.e. cells), resolving kinetics, conformational changes ...



* Adato et al., *Nat. Comm.* **4**, 2154 (2013)

spectral overlap

• Complex sample analysis require chemometrics





Surface-enhanced IR absorbance spectroscopy (SEIRA)




IR absorbance \propto Oscillator strength of molecular bond x | Electric field | 2

Optical nanoantennas can be used to boost IR absorbance






Signal



Surface-enhanced IR absorbance spectroscopy (SEIRA)





- IR absorbance \propto Oscillator strength of molecular bond x | Electric field | 2
- Optical nanoantennas can be used to boost IR absorbance.




- Mid-IR spectroscopy is broadband & provides rich information
- If a single/narrow resonance is used, collected spectral information will be limited


Life is formed by Biomolecules from Four Major Classes






Real-time SEIRA measurements in aqueous media - nanoplasmonic metasurfaces










Nature Communications 2013, patent granted Nano Letters 2016 Nature Communications 2018 **Advanced Materials, 2021**

Heterogeneous bioassay in solution







Multi-resonant metasurfaces tuned for lipid, peptide, nucleic acids, polysaccharide




Nature Communications 2013, patent granted Nano Letters 2016 Nature Communications 2018 **Advanced Materials, 2021**

Heterogeneous bioassay in solution






Real-time & label-free observation of vesicle-peptite interactions & cargo release





Reminder: Mid-infrared sensing for biology




non invasive access to molecular vibrations


- signal intrinsic to molecular structure
- Resolves secondary structure from amide I band shape







Protein structure













Misfolding disrupts protein function & lead to diseases




Protein misfolding – a common mechanism in neuordegenerative disorders






Aggregation pathway of alpha-synuclein protein



ImmunoSEIRA for protein structure analysis


Sensing scheme




- Real time label-free optofluidic setup
- Immunoassay captures target proteins : α- syn, with highest specificity irrespective of their conformation using antibodies
- From the Amide I absorption bands (1600-1700 cm⁻¹) of the α- syn species, their distinct secondary structures are retrieved

Optofluidic configuration


2D microarray SEIRA design




Small sample volume microfluidics





Optofluidic SEIRA signal




Conformational profiling of different alpha-synuclein protein forms











But, how can we quantitatively separate the pathological aggregates-oligomers and fibrils in a mixed sample?

Al-Boosted Immuno-SEIRA for quantitative diffrentiation of aggregate mixture




Multiplexed Detection of Structural Biomarkers with Immuno-SEIRA Microarray





Confidential



Wavenumber (cm-1)

110

130


150

170


Time (minutes)

190


210





Signal obtained on SYN-211 antibody functionalized unit Integrated Absorbance (mODcm⁻¹) **Butter wash** 200 100 150 Time (minutes) **Blocking wash** Absorbance (mOD) Tau & ASyn fibrils injection 1500 1600 1700



How to use high-Q factor Mid-IR antennas?







PRL 101, 157403 (2008)


J. Am. Chem. Soc. 135, 3688–3695 (2013)

Nat. Commun. 9, 2160 (2018)


State-of-the art in SEIRAS - plasmonic antennas




limited by broad resonances




- What can we do differently with sharper resonances?
- How can we harness them for new functionalities?




Mid-IR imaging biosensor with high-Q dielectric metasurfaces





Nanofabrication of pixelated dielectric metasurfaces



100 metapixels, total size: 2.25 mm²

Fabrication method:

- hydrogenated amorphous silicon layer deposition via PECVD
- MgF₂ substrate
- electron-beam lithography
- fluorine-based dry plasma etching

linear scaling of resonator size, S


A = 1.96 um

B = 0.96 um


 $P_{x} = 3.92 \text{ um}$

 $P_{y} = 2.26 \text{ um}$

Fixed height = 0.7 um Fixed tilt angle θ = 20°




Broadband & Imaging-Based Spectral Detection : Molecular Barcodes




Tittl et al., Science **360**, (2018)


Chemical-specific barcodes: image based identification


Biosensing & diagnostics






Materials science





Environmental monitoring





Barcode images are unique & different




Absorption spectrum of each sample is distinct

Tittl et al., Science 360, (2018)

Broadband Mid-IR Metasurfaces for SEIRA Spectroscopy & BioSensing


Pixelated Metasurfaces:

- pixel size >100µm
- discrete spectral sampling
- → large area required




Gradient Metasurface Concept:

- resonances tolerate perturbation
- scale every column slightly
- → continues spectral coverage
- → much less area required




Broadband Gradient Mid-IR Metasurfaces




Broadband Gradient Metasurface for SEIRA Spectroscopy & BioSensing

Single-Analyte: Detection

Metasurface is coated with a polymer, such as PMMA, PS or PVA






Broadband Gradient Metasurface for SEIRA Spectroscopy & BioSensing


Single-Analyte: Detection

Metasurface is coated with a polymer, such as PMMA, PS or PVA





Multi-Analyte: Decomposition


Metasurface is coated with a polymer mixture




Flexibilities Offered by Gradient Mid-IR Metasurfaces













High Harmanic Generation: Nonlinear Optics with Gradient Broadband Metasurfaces

In collaboration with Yuri Kivshar's Group at ANU











Jangid & Richter et al., Advanced Materials (2023)

Wafer-scale fabrication of mid-IR metasurfaces







- 90% transmission over an ultra-broadband spectrum [from 2 to 20 μm]
- Ultra-low effective refractive index down to 1.04





Wafer-scale fabrication of mid-IR metasurfaces











A. Leitis et.al., Advanced Materials, 2102232 (2021)

BioNanoPhotonic Systems Labotatory



Ongoing Research Directions








Live Cell & Single-Cell Technologies

- Cell screening & profiling, e.g. for immunoengineering applications
- Real-time & label-free secretion analysis from single cells
- Cell-cell interactions & incorporation of biomimetic surfaces








ACS Sensors 2018 Small 2018 Science Advance 2021 Biosensors and Bioelectronics 2022 Nature Biomedical Engineering 2023 ...

Surface Enhanced IR Spectroscopy

- Chemical-specific & fingerprint detection
- Monitor conformational/structural changes
- Real-time measurements in aqueous solution



Since PNAS 2009 Recent highlights:

Nature Communications 2013 → patent granted Science 2015

Nature Communications 2018

Science 2018


Science Advance 2019

Advanced Materials 2021 & 2022

PNAS 2023 ...

Point-of-Care Diagnostics

- Low-cost, portable, rapid & multiplexed microarray technologies.
- New device schemes including lens-free microscopy, hyperspectral imaging ...



Since PNAS 2011 ... Recent highlights

Light Science and Technology 2018

ACS Photonics 2015 → patent granted

ACS Nano 2018 → patent filed


Nature Photonics 2019


Small 2020

Nature Communications 2021


Importance of heterogeneity

Cells present heterogeneity in various characteristics






Homogenous population




Heterogeneous population




Highly heterogeneous population


Protein secretion from cells


- Closely related to physiology and disease progression
- A precisely regulated dynamic process
 - External stimulation
 - Gene expression
 - Exocytosis
- Involved in biological functions, interactions & communications
- A variety of secreted substances:
 - Cytokines
 - Antibodies
 - Extracellular vesicles ...



Single-Cell Microarray for "Spatio-temporal" Secretion Analysis









(July Issue 2023)

Nanoplasmonic single-cell microwell arrays









Integrated biochip:

- Maintains live cells



PDMS microwells:

- Compartmentalize single cells





Plasmonic nanoholes


- Enables label-free detection


Label-free Detection Principle - NanoPlasmonics

Plasmonic Gold Chip with Nanohole Arrays










Confined near-fields within <50 nm from surface enables:


- → Strong analyte-light overlap
- → High sensitivity


High-throughput spatiotemporal monitoring of single-cell secretions




Deterministic single-cell loading in the microarray









Specific detection of target analyte – surface functionalization








Wafer-scale manufacturing of nanoplasmonic sensors with DUVL process




- Substrate is typically glass or silicon
- This example is based on glass:
- Coated with 100 nm thick gold film
- ~50 1x1 cm² nanoplasmonic chips



Simultaneous monitoring of cell secretion and morphology



Spatio-temporal visualization of extracellular secretion from single cells




Max ≡ **Maximum amount of secretion**

Slope **≡ Secretion rate**


Max ≡ Maximum areal coverage

Slope **≡ Adsorption rate**

Statistical analysis of a population with single-cell microarray




Protein transport inhibited (PTI) hybridoma








Cell secretion during mitosis

Non-destructive monitoring of secretion during cell division










Downregulation of the ER-Golgi pathway before mitosis

Cell content release during cell death






- cell shrinkage
- plasma membrane blebbing
- formation of apoptotic bodies



necroptosis




- increase in cell volume
- loss of plasma membrane integrity
- leakage of cellular contents


Different forms of cell deaths:

- Programmed/regulated vs unregulated
- Apoptosis vs necrosis

Cell content release during apoptosis and necroptosis




Apoptosis is immunologically "silent"




Necroptosis

Bursting behavior of necroptosis could mediate:

- "long-distance" cell communications
- "Higher immunogenicity"




EPFL Photonic biosensors & bioimaging platforms for cell screening




- Cell-cell interactions for cytotoxicity analysis
- Spatio-temporal secretion analysis


Cytotoxic T cells in anti-tumor immunity



Cytotoxic T cells



High-throughput biochip for real-time analysis of cell-cell interactions






Single Cell Microarrays For Monitoring of Cancer Cell Killings



Time-lapse imaging with 4 simultaneous channels (1 bright field & 3 fluorescence)







Absence of "cell death"

Specific "death"
Caspase 3/7
pathway


Bright field image




T-cells



Tumor cells



"Cytotoxicity" Effect of T Cells on Melanoma Cancer Cells




Technology generates a large set of spatially and temporally resolved data

→ Images of each cell from multiple channels (bright-field & fluorescent channels) over a long time period (i.e. 24 hours with 15 min intervals)



Deep neural network:

Rapid and quantitative analysis of single cell killing events at high-throughput







BIOS: BioNanoPhotonic Systems Labotatory

Mission

Develop science and technology for powerful optical bioanalytical devices that can have impact on areas ranging from basic research in life sciences to real-world applications such as disease diagnostics.




Microfluidics



Nanofabrication








Data Science

Acknowledgements





Scientific Collaborators:

- Hilal Lashuel, EPFL
- Camilla Jandus, University of Geneva
- Alexander Mathis, EPFL
- Matthias Lutholf, EPFL & Roche
- George Coukos, Agora
- Yuri Kivshar, Australian National University
- Javier Garcia de Abajo, ICFO, Spain
- Sang-Hyun Oh, Univ. Minnesota
- ...and more

Open MS thesis projects!!

Check our website & send your CV bios.epfl.ch hatice.altug@epfl.ch

Open post-doc and PhD positions!!

Send your CV bios.epfl.ch hatice.altug@epfl.ch



Acknowledgements



- Swiss State Secretariat for Education, Research and Innovation (SERI)
 - under contract number 22.00081 for **MIRAQLS** project (#101070700) linked to European Union Horizon Europe Research and Innovation Programme
 - under contract number 22.00018 for **TWISTEDNANO** project (#101046424) linked to European Union Horizon Europe Research and Innovation Programme
- European Union's Horizon 2020 Research and Innovation Programme
 -under the Marie Skłodowska-Curie grant agreement No. 955623 (CONSENSE Project)
- Swiss National Science Foundation Korean-Swiss Science & Tech Programme
- Swiss National Science Foundation Sinergia Grant
- EPFL & ZEISS IDEAS Program



THANKYOUI

Prof. Hatice ALTUG

BioNanoPhotonic Systems Laboratory

École Polytechnique Fédérale de Lausanne, Switzerland