






Applications Lesson 12

MSE 304


Nako Nakatsuka


Nako.Nakatsuka@epfl.ch




Recap of Last Week's Class on Applications

- Surface Chemistry alkanethiols, alkylsilanes, polymers for nonspecific binding
- Nanoscience interactions of nanoparticles with light (plasmon resonance)
- AFM and STM to Visualize Surfaces












What We Cover in Today's Class

- Plasmonic nanoparticles short recap on plasmon resonance used for rapid tests!
- What is a Biosensor molecular recognition and performance parameters
- History of Biosensors Clark electrode and recap of electrochemistry
- State-of-the-art Biosensor: Glucose Sensor recap of electrolysis
- Remaining challenges of Biosensors now let's take a look at the COVID-19 test




Influencing the Sea of Electrons on Surfaces Using Light

The Free Electron Sea



Light as an oscillating electromagnetic field




Copyright 1998 by John Wiley and Sons, Inc. All rights reserved.

Electrons in metal are affected by the electromagnetic field from light incident to surface




Effect of Plasmon Resonance



Incident light excites electrons on metal surface

Oscillating field of light interacts with free electrons





Electron cloud is perturbed

Electrons begin to oscillate collectively as a polarized cloud (plasmon)



Effect of Plasmon Resonance







Electron cloud is perturbed

Electrons begin to oscillate collectively as a polarized cloud (plasmon)



Effect of Plasmon Resonance





At resonance frequency, incident light energy transferred to oscillating electron cloud (plasmon excitation)

Only specific wavelengths of light are reflected back while others are absorbed or scattered

Relative contributions of absorption and scattering determined by size and shape of plasmonic particles



Going Back Many Centuries (4th Century): Lycurgus Cup

The glass contains gold and silver nanoparticles 10–100 nm in size

Resonance scatters green wavelengths



Lit from front (reflection)




(transmission)

Nanoparticles absorb shorter wavelengths (blue and green) and remaining wavelengths pass through and is transmitted




COVID Rapid Tests – Use of Gold Nanoparticles




First, we need to understand how biosensors work!



Biosensors For Human Health




Molecular Recognition Elements/Biological Receptors






How do Bioreceptors Interact with their Analytes?

Expectation




Molecular recognition of target by bioreceptor






Biosensor: Transduces Molecular Binding Event to a Signal





How Do We Smell Things?





The World's Best "Biosensor"

Olfactory epithelium

Back of nasal passage that houses sensory receptors

Humans 5-10 cm²

Dogs 60-170 cm²

Olfactory bulb

Brain region responsible for processing smells

Dogs have **40x** larger olfactory bulbs than humans relative to total brain size



Olfactory receptors

Back of nasal passage that houses sensory receptors

Humans 5-6 million

400 different types of receptor proteins


Dogs 220-300 million

800-1200 different types pf receptor proteins


Dogs are estimated to be 10,000 – 100,000 x more sensitive to certain odors than humans




Mimicking Nature: Deep-Learning Inspired Electronic Nose




How to Read Out the Signal of a Biosensor








How to Read Out the Signal of a Biosensor









Concentration – What Do these Values Mean?

 $1 M = 1 mol/L = 6 x 10^{23} molecules/L$

Note: Avogadro's number: 6 x 10²³ molecules/mol


1 M \sim 1 x 10¹⁵ molecules/100 μ m³









How to Read Out the Signal of a Biosensor








How to Calculate the Equilibrium Binding Affinity (K_d)




$$\phi = [AR]$$
 Target-bound receptors $[R_0]$ Total number of receptors






How to Calculate the Sensitivity and Limit of Detection (LOD)


 $m LOD = rac{3\sigma_{S_0}}{rac{dS}{dS}}$ 3 x standard deviation of zero analyte signal









The Detection Range is Dependent on the K_d of the Bioreceptor







Different Signal Transduction Mechanisms of Biosensors






History of Biosensors



Leland Clark The "Father of Biosensors"




Electrochemistry: Relationship between Chemical Reactions + Electricity

Electricity: movement of electrons


 $e^- \longrightarrow e^- \longrightarrow e^-$

Certain chemical reactions -> create electricity




Oxidation & Reduction Reactions

Spontaneous reaction




Electricity -> make certain chemical reactions happen




Non spontaneous reaction



Galvanic/Voltaic Cells for Electrochemistry



Standard Reduction Potential of Different Metals




Stronger pull for electrons

Under standard conditions:

- 25 °C (298 K)
- 1 M concentration for reactants/products
- Relative to the standard hydrogen electrode (0.00 V)


$$\mathrm{H_2} \leftrightarrow 2\mathrm{H}^+ + 2e^-$$



Cu will pull electrons from Zn



Galvanic/Voltaic Cells for Electrochemistry




Clark Electrode Measures Ambient Oxygen Partial Pressure

Clark Electrode was originally designed to measure oxygen concentration in blood



Leland Clark The "Father of Biosensors"






Clark Electrode Measures Ambient Oxygen Partial Pressure



Leland Clark The "Father of Biosensors"






Clark Electrode Measures Ambient Oxygen Partial Pressure

Redox processes are highly favorable so spontaneous process


Electron flow between electrodes generates current directly proportional to the partial pressure of oxygen in the sample








Glucose Biosensor – Transformative for Human Health










Todaro *et al.*, *Front. Chem., 10*, **2022**

Glucose Sensor – Incorporating an Enzyme (Biological Catalyst)

$$\begin{array}{c} \text{OH} \\ \text{HO} \\ \text{OH} \\ \text{D-glucose} \\ \text{D-gluconolactone} \end{array} + \begin{array}{c} \text{Glucose} \\ \text{oxidase} \\ \text{D-gluconolactone} \\ \end{array} + \begin{array}{c} \text{OH} \\ \text{HO} \\ \text{OH} \\ \text{$$


Catalyst: substance that accelerates the reaction rate without being consumed in the process

Catalysis: affects kinetics and not thermodynamics (it makes kinetics faster by decreasing activation energy)







Intermolecular Interactions Between Enzyme and Substrate




Active site of glucose oxidase



Bauer, *Biomolecules*, 12, 3, 2022











How Can Enzymes Catalyze Reactions?




Enzymes lower the activation energy of chemical reactions, allowing reactions to occur faster at the temperatures and conditions found in biological systems.







Glucose Sensor – Detection of Hydrogen Peroxide as Readout



Monitor decrease in O₂ as proxy for glucose

Lack of **sensitivity** – small signal change since O₂ levels high in biofluids

$$S = rac{\Delta ext{Signal}}{\Delta ext{Concentration}}$$

Lack of **specificity** – O₂ sensors influenced by changes unrelated to glucose (respiratory, temperature)


Monitor increase in H₂O₂ as proxy for glucose

Higher sensitivity – H₂O₂ not abundant in the sample prior to enzymatic reaction

Higher specificity – H₂O₂ is direct product of enzymatic reaction



Oxidation of H₂O₂ on the Surface of Platinum Electrodes



Half Reaction						potential
	F ₂	+	2e ⁻	\rightleftharpoons	2F ⁻	+2.87 V
	Pb ⁴⁺	+	2e ⁻	\Longrightarrow	Pb ²⁺	+1.67 V
	CI ₂	+	2e ⁻	\Longrightarrow	2Cl ⁻	+1.36 V
	Ag⁺	+	1e⁻	\Longrightarrow	Ag	+0.80 V
	Fe ³⁺	+	1e⁻	\Longrightarrow	Fe ²⁺	+0.77 V
	Cu ²⁺	+	2e ⁻	\Longrightarrow	Cu	+0.34 V
	2H⁺	+	2 e ⁻	\rightleftharpoons	H ₂	0.00 V


Standard reduction potential for H_2O_2 : $E^0 = +0.695$ V

So $E_{\text{oxidation}} = -0.695 \text{ V}$


The oxidation process is non-spontaneous!




Electrolysis - Driving Chemical Reactions with Electricity



Water splitting



Oxidation




Non-spontaneous reaction


		ı	Half F	React	ion		potential
		F ₂	+	2e ⁻	\Longrightarrow	2F-	+2.87 V
	1	Pb⁴⁺	+	2e ⁻	\Longrightarrow	Pb ²⁺	+1.67 V
		O ₂	+	2e ⁻	+ 4H+	\rightleftharpoons	2H ₂ O +1.23 V
		Ag⁺	+	1e⁻	\Longrightarrow	Ag	+0.80 V
Stronger pu		Fe ³⁺	+	1e⁻	\Longrightarrow	Fe ²⁺	+0.77 V
for electror	1S	Cu ²⁺	+	2e ⁻	\Longrightarrow	Cu	+0.34 V
		2H+	+	2e ⁻	\Rightarrow	H ₂	0.00 V

We can force this reaction to happen using electricity



Electrolytic Cell - Driving Chemical Reactions with Electricity






Applying a potential drives non-spontaneous redox reaction


Operating potential applied to overcome activation energy barrier to drive reaction

CHEMINA


Glucose Sensor Requires an Operating Potential



Success of the Glucose Biosensor - Robust System



Highly specific and robust enzyme



Detection window for H₂O₂ doesn't overlap with other species in blood



 H_2O_2




LEVEL	mg/dl	mmol/L	RISK	SUGGESTED ACTION
DANGER - HIGH	315+	17.4	VERY HIGH	MEDICAL ATTENTION
HIGH	280	15.6	HIGH	MEDICAL ATTENTION
HIGH	250	13.7	HIGH	MEDICAL ATTENTION
HIGH	215	11	HIGH	MEDICAL ATTENTION
BORDERLINE	180	10	MEDIUM	CONSULT DOCTOR
BORDERLINE	150	8.2	MEDIUM	CONSULT DOCTOR
BORDERLINE	120	7	MEDIUM	CONSULT DOCTOR
NORMAL	108	6	NO RISK	NO ACTION NEEDED
NORMAL	72	4	NO RISK	NO ACTION NEEDED
LOW	70	3.9	MEDIUM	CONSULT DOCTOR
DANGER - LOW	50	2.8	HIGH	MEDICAL ATTENTION

Clear values associated with healthy and disease states (diabetes)




Success of the Glucose Biosensor - Concentration Range






1. Nonspecific Binding (NSB) in Complex Environments

Expectation (phosphate buffered saline)




Reality (biofluids such as blood)





Most Existing Diagnostics are Limited by Nonspecific Binding


Non-specific binding induced noise at low concentrations limits diagnostic biosensors






2. Very Few Number of Molecules Within Sample Size

Stochastic process – probability of the biosensor seeing such a low concentration of target diffusing to the surface?







Few Molecules in Large Volume - Diffusion to the Sensor Surface

Stochastic process – probability of the biosensor seeing such a low concentration of target diffusing to the surface?


100 μ m 1 μ M ~1 billion 1 pM ~1000 1 fM ~1 Single molecule of salt in several Olympic-sized swimming pools!

Making it stick: convection, reaction and diffusion in surface-based biosensors

Todd M Squires ☑, Robert J Messinger & Scott R Manalis


based on diffusion (D)

Nature Biotechnology 26, 417–426 (2008) Cite this article





3. Challenges of Bioreceptor Specificity in Addition to NSB




Bioreceptors Can be Cross-Reactive to Nonspecific Molecules

Reproducibility crisis: Blame it on the antib

Monya Baker

Nature **521**, 274–276 (2015) | Cite this article



Cross-reactivity: binding to other non-targets

Batch-to-batch variability


Sensitivity to environment






4. Challenges of Continuous Monitoring








Continuous Measurements: Reversible Binding to Bioreceptors



We need to engineer/redesign bioreceptors to overcome this mismatch



5. Tissue Response to Implanted Sensors Prevents Long-Term Use






Size Matters: Miniaturization of Sensors

Microglia activation hours after implantation

Carbon fiber electrodes (20-60 µm²)




Microdialysis probes (30,000-70,000 μm²)


Planar silicon electrodes (3000-15,000 μm²)



Let's Go Back to The COVID-19 Test: Antibody Sandwich Assays



Visualize binding of gold nanoparticles (AuNP) by eye




Nucleocapsid protein in high abundance and stable in nasal swabs






Antibody Sandwich Assays are Commonly Used in Biosensors




Visualize binding of gold nanoparticles (AuNP) by eye



Need a detector (expensive) to visualize fluorescence




Challenges of Antibody Sandwich Assays



Steric hindrance at the surface

Density of bioreceptors need to be optimized

False negatives (not enough AuNPs bound)




Nonspecific sticking of the gold nanoparticles

To the surface, to primary antibodies, etc.


False positives (too many AuNPs bound)



The COVID-19 Test Uses and Antibody Sandwich Assay




Steric hindrance at the surface



Nonspecific sticking to gold nanoparticles




Lateral Flow Assay Optimized to Detect COVID-19



Yes/No test – not quantitative!



Next Week: DNA Aptamers for Biosensing on Semiconductors





Summary of Today's Class

- Plasmonic nanoparticles short recap on plasmon resonance used for rapid tests!
- What is a Biosensor molecular recognition and performance parameters
- History of Biosensors Clark electrode and recap of electrochemistry
- State-of-the-art Biosensor: Glucose Sensor recap of electrolysis
- Remaining challenges of Biosensors
- How the COVID-19 rapid test works

Exercise Session: Real-world applications of biosensors for human health (glucose and neurotransmitters)

