

- Installation of Dr. Probe¹
 - Graphical user interface (GUI) for Windows: [installer download](#) and [system requirements](#)
 - Unzip and run the .exe

- No need of license activation for academics

- Initial step with the GUI
 - Start Dr. Probe: start menu / ERCtools / Dr. Probe
 - Set the preferences: Setup / Preferences
 - select your GPU device
 - enter the number of CPU, according to the guidelines (~half of the CPU number of your computer that you can read in the “Status information” window)
 - Save these parameters: File / Exit

- Documentation
 - [Website](#)
 - [Paper](#)