

- Installation of Dr. Probe (GUI for Windows)
 - Download the Windows installer: <https://pgi-jcns.fz-juelich.de/barthel/drprobe/drprobe-download.html>
 - Unzip and run the .exe
- Activate license → free for academics
 - Start Dr. Probe GUI (start menu / ERCtools / Dr. Probe)
 - Enter a username of your choice
 - Enter the following registration string: *2023-summer-break-trial*
 - Press register (200 days trial)
 - For a permanent license, send your request string to drprobetool@gmail.com to get your registration string
- Initial step in the GUI
 - Setup / Preferences: select the GPU and enter the number of CPU to be used
- Documentation
 - Barthel, *Ultramicroscopy* **193**, 2018: <https://doi.org/10.1016/j.ultramic.2018.06.003>
 - Website: <https://pgi-jcns.fz-juelich.de/barthel/drprobe/index.html>